


Custom Peptide Service

Any Sequence, Any Length Any Complexity



GenScript's proprietary FlexPeptide™ platform provides custom peptides of up to 200 residues in length with a capacity of 10,000 peptides per month. With a scale ranging from milligrams to kilograms, comprehensive labeling and modification options, and rapid high-throughput peptide library synthesis, our platform is an ideal choice for peptide synthesis, peptide library, and peptide array services. GenScript also offers large-scale cosmetic peptide and cGMP-grade peptide production services to accommodate both product manufacture and drug development, respectively.

◆ Peptide Synthesis

**Any sequence, any length, any complexity;
Comprehensive modifications; 95% success rate**

- TQM system and ISO9001:208 certified
- Online quotation and order
- Free peptide antigen design

◆ Peptide Array

**For small quantities of 48 peptides or more for
protein-protein interaction and enzyme substrate
screening**

◆ Cosmetic Peptide Synthesis

**Large-scale, economically priced, cosmetic
peptide manufacturing**

◆ Click Peptide Synthesis

**Special chemistry for soluble β -amyloid
synthesis**

◆ Peptide Library

**For orders of 24 peptides or more; up to 10,000
peptides per month for vaccine development,
epitope mapping, and drug discovery**

- Cost-efficient pricing
- Epitope mapping downstream service

◆ Large-Scale Peptide Synthesis

**Multi-gram level custom peptide manufacturing
for any application**

◆ cGMP Peptide Synthesis

**Fully documented, confidential service for the
development of peptide drugs and vaccines**

◆ Peptoid Synthesis

Custom peptidomimetics for drug discovery

Comprehensive Modifications

**Including Amidation and Acetylation • Biotinylation • Phosphorylation • Methylation • FITC labeling •
PEGylation • KLH, BSA, and OVA conjugation**

www.genscript.com

Total Peptide Quality Management

The most stringent QC/QA in the industry

Total Quality Management

With HPLC and MS analyses after every peptide preparation step, GenScript's Total Quality Management (TQM) platform ensures that your custom peptide is delivered with the most accurate purity, and highest quality possible.



AccuPep+ QC




AccuPep+ QC is a collection of additional quality control services, such as solubility testing and TFA counter-ion removal, that can ensure the future reliability and reproducibility of your custom peptides in cellular assays, kinetic studies, and more.

ArgonShield™




Every custom peptide is protected from oxidation and deliquescence by our ArgonShield™ packing service, maintaining the integrity of your custom peptide from our lab to yours.

Peptide Technical Resources




Bioinformatics
Tools



Protocols



Frequently Asked
Questions



Peer-reviewed
Publications

www.genscript.com

To request a quote online visit www.genscript.com and click: **Services** ► **Peptide Services** ►
To request a quote by email send your sequence, purity, and quantity to peptide@genscript.com

