

Human Recombinant 4-1BB Stable Cell Line
Cat. No. M00609**Version 04282015****I. INTRODUCTION**

Catalog Number: M00609

Cell Line Name: GS-H2/4-1BB

Gene Synonyms: 4-1BB; CD137; ILA; TNFRSF9

Expressed Gene: Codon Optimized from NM_001561.5; no expressed tags

Host Cell: GS-H2

Quantity: Two vials of frozen cells (1×10^6 per vial)

Stability: 20 passages

Application: *in vitro* functional assay

Freeze Medium: 95% complete growth medium, 5% DMSO

Complete Growth Medium: MEM, 10% FBS

Culture Medium: MEM, 10% FBS, 200 $\mu\text{g/ml}$ Hygromycin B, 2 $\mu\text{g/ml}$ Puromycin

Mycoplasma Status: Negative

Functional Performance: For 4-1BBL (ligand of 4-1BB), Signal / Background (S/B) > 3

Storage: Liquid nitrogen immediately upon receipt

II. BACKGROUND

4-1BB is a member of the tumor necrosis factor (TNF) receptor family. Its alternative names are tumor necrosis factor receptor superfamily member 9 (TNFRSF9), CD137 and induced by lymphocyte activation (ILA). It is currently of interest to immunologists as a co-stimulatory immune checkpoint molecule.

4-1BB can be expressed by activated T cells, but to a larger extent on CD8 than on CD4 T cells. In addition, CD137 expression is found on dendritic cells, follicular dendritic cells, natural killer cells, granulocytes and cells of blood vessel walls at sites of inflammation.

The best characterized activity of 4-1BB is its costimulatory activity for activated T cells. Crosslinking of 4-1BB enhances T cell proliferation, IL-2 secretion, survival and cytolytic activity. Further, it can enhance immune activity to eliminate tumors in mice.

III. REPRESENTATIVE DATA

- Protein Expression Validation




Figure 1. Flow cytometry analysis of 4-1BB protein expression in GS-H2/4-1BB cells. Red: GS-H2, Blue: GS-H2/4-1BB.

- Validation by *in vitro* Functional Assay




Figure 2. Functional analysis of 4-1BBL on GS-H2/4-1BB cells by measuring IL-8 release.

- Stability Validation

	G1	G5	G15	G15	G20
EC50 (µg/ml)	0.4793	0.2525	0.2731	0.3374	0.3263
S/B	5.14	4.37	3.40	4.28	4.80

IV. THAWING AND SUBCULTURING

Thawing Protocol

1. Remove the vial from liquid nitrogen tank and thaw cells quickly in a 37°C water-bath.
2. Just before the cells are completely thawed, decontaminate the outside of the vial with 70% ethanol and transfer the cells to a 15 ml centrifuge tube containing 9 ml of complete growth medium.
3. Pellet cells by centrifugation at 200 x g for 5 min, and remove the medium.
4. Resuspend the cells in complete growth medium.
5. Transfer the cell suspension to a 10 cm dish with 10 ml of complete growth medium.
6. Grow the cells in incubator with 37°C, 5 % CO₂.
7. Add antibiotic in the following day.

Sub-culturing Protocol

1. Centrifuge the cells at 200 x g for 5min, and remove the medium.
2. Resuspend the cells in culture medium and add the cells suspension to new culture dish.
3. Grow the cells in incubator with 37°C, 5 % CO₂.

Subcultivation Ratio: A subcultivation ratio of 1:4 to 1:8 is recommended

Medium Renewal: Every 2 to 3 days

V. REFERENCES

1. Kwon BS, Weissman SM (1989). "cDNA sequences of two inducible T-cell genes". *Proc. Natl. Acad. Sci. U.S.A.* 86 (6): 1963–7.
2. Schwarz H, Tuckwell J, Lotz M (1993). "A receptor induced by lymphocyte activation (ILA): a new member of the human nerve-growth-factor/tumor-necrosis-factor receptor family". *Gene.* 134 (2): 295–8.
3. Sica G, Chen L (2000). "Biochemical and immunological characteristics of 4-1BB (CD137) receptor and ligand and potential applications in cancer therapy". *Arch. Immunol. Ther. Exp. (Warsz.)*. 47 (5): 275–9.
4. Schwarz H (2005). "Biological activities of reverse signal transduction through CD137 ligand". *J. Leukoc. Biol.* 77 (3): 281–6

GenScript USA Inc,
860 Centennial Ave.
Piscataway, NJ 08854
Toll-Free: 1-877-436-7274
Tel: 1-732-885-9188, Fax: 1-732-210-0262
Email: product@genscript.com
Web: <http://www.genscript.com>
For Research Use Only.

Limited Use License Agreement

This is a legal agreement between you (Licensee) and GenScript USA Inc. governing use of GenScript's stable cell line products and protocols provided to licensee. By purchasing and using the stable cell line, the buyer agrees to comply with the following terms and conditions of this label license and recognizes and agrees to such restrictions:

- 1) The products are not transferable and will be used at the site where they were purchased. Transfer to another site owned by buyer will be permitted only upon written request by buyer followed by subsequent written approval by GenScript.
- 2) The purchaser cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party.
- 3) The products sold by GenScript are for laboratory and animal research purposes only. The products are not to be used on humans, for consumption, or for any unlawful uses.

GenScript USA Inc. will not assert against the buyer a claim of infringement of patents owned or controlled by GenScript USA Inc. and claiming this product based upon the manufacture, use or sale of a clinical diagnostic, therapeutic and vaccine, or prophylactic product developed in research by the buyer in which this product or its components has been employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on the use of this product for other purposes, contact Marketing Department, GenScript USA Inc., GenScript 860 Centennial Ave Piscataway, NJ 08854. Phone: 1-732-885-9188. Fax: 1-732-210-0262. Email: marketing@genscript.com.