



MCS-1: 4525-4719

Polyhedrin promoter

--- TTTTACTGTT TTCGTAACAG TTTTGTAAAT AAAAAAACCTATA AATATTCC GGATTATTCA TACCGTCCCA

CCATCGGGCG CGGATCCCGG TCCGAAGCGC GCGGAATTCA AAGGCCTACG TCGACGAGCT CACTAGTC

BamH I *EcoR I* *Stu I* *Sal I* *Sac I* *Spe I*

GC GGCCGCTTTC GAATCTAGAG CCTGCAGTCT CGACAAGC TTGTCGAGAA GTACTAG---

Not I *BstB I* **Stop** *Pst I* *Hind III* **Stop**

MCS-2: 4390-4251

p10 promoter

---ATATTAATTA AAATACTATA CTGTAAATTA CATTTTATT ACAATCACTC GACGAAGACT TGATCACCCG GG

Sma I

ATCTCGAG CCATGGTGCT AGCAGCTGAT GCATAGCATG CGGTACCGGG AGATGGGGGA GGCTAACTGA---

Xho I *Nco I* *Nhe I* *Sph I* *Kpn I* **Stop** **Stop**