



**Polyhedrin promoter**

--- TTTTACTGTT TTCGTAACAG TTTTGTAAT AAAAAAACCTATA AATATTCC GGATTATTCA TACCGTCCCA  
 CCATCGGGCG CGGATCCCGG TCCGAAGCGC GCGGAATTCA AAGGCCTACG TCGACGAGCT CACTAGTC  
 GC GGCCGCTTTC GAATCTAGAG CCTGCAGTCT CGAGGCATGC GGTACCAAGC TTGTCTGAGAA GTACTAG---

Restriction enzymes for MCS:

- BamHI*
- EcoRI*
- SalI*
- SacI*
- SpeI*
- NotI*
- BstBI*
- XbaI*
- PstI*
- XhoI*
- KpnI*
- HindIII*