


Summary:
 Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient : A B
 0.1min 75 25
 25 min 50 50
 25.01min 0 100
 30 min 0 100
 Flow rate : 1.0ml/min

Cat. No.: RP11056
 Lot. No.: PX00228502



RESULTS

| Peak No. | Peak ID | Ret Time | Height | Area | Conc |
|----------|---------|----------|------------|-------------|----------|
| 1 | | 9.617 | 1031.749 | 30247.111 | 0.4825 |
| 2 | | 10.548 | 352.234 | 4307.102 | 0.0687 |
| 3 | | 11.220 | 171854.688 | 6234412.000 | 99.4488 |
| TOTAL | | | 173238.670 | 6268966.213 | 100.0000 |