


Mass Spectrum



Sample Information

|                         |            |                 |                 |
|-------------------------|------------|-----------------|-----------------|
| Month-Day Processed :   | 12/09/25   | Equipment       | :ZJ22010150     |
| Time Processed :        | 14:25:34   | Interface Bias  | : +4.5 kV       |
| Injection Volume :      | 0.2        | Drying Gas Flow | :5 L/min        |
| Sample Name :           | Flagelin22 | T.Flow          | :0.2 ml/min     |
| Cat No. : RP19986       |            | B.conc          | :50%H2O/50%MeOH |
| Lot No. : E0283025K     |            |                 |                 |
| Theoretical MW: 2272.50 |            |                 |                 |
| Observed MW : 2272.5    |            |                 |                 |