

Sample Name:C-Peptide (3-33), human
Cat No. : RP11116
Lot No. : E0597726C
Time Processed :16:42:56
Month-Day-Year Processed :03/19/2026


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ21010376

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	15.950	9270971	1557964	97.093
2	16.150	190873	28892	1.999
3	16.258	30645	9677	0.321
4	16.383	56069	5779	0.587
Total		9548559	1602312	100.000