

Sample Name:C-Peptide (3-33), human

Cat No. : RP11116

Lot No. : E0980225J

Time Processed :8:20:17

Month-Day-Year Processed :10/15/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	11.847	10698	1433	0.146
2	12.119	35407	1502	0.483
3	15.361	7226393	1025169	98.478
4	16.220	65577	3374	0.894
Total		7338075	1031478	100.000