

Sample Name :KRAS G12D Peptide (VVVGADGVGK)

Cat No. : RP30227

Lot No. : E1265425E

Time Processed :11:17:56

Month-Day-Year Processed :05/20/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ23010113

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.273	5000	1124	0.833
2	8.465	2311	548	0.385
3	8.620	570302	110727	95.005
4	8.785	6493	1391	1.082
5	8.997	16179	3648	2.695
Total		600285	117439	100.000