

Sample Name :C-Peptide (3-33), human
Cat No. :RP11116
Lot No. :E1804225H
Time Processed :15:26:53
Month-Day-Year Processed :09/02/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ19010142

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	13.242	23052	3399	0.647
2	13.400	3468374	531573	97.295
3	13.550	36972	11479	1.037
4	14.492	34216	6124	0.960
5	15.933	2204	321	0.062
Total		3564818	552896	100.000