


Mass Spectrum



Sample Information		Interface	:ESI	Equipment	:ZJ22010150
Month-Day Processed :	09/02/25	Nebulizing Gas Flow	:1.5L/min	Interface Bias	: +4.5 kV
Time Processed :	1:36:19 PM	CDL Temp	:250	Drying Gas Flow	:5 L/min
Injection Volume :	0.4	Block Temp	:200	T.Flow	:0.2 ml/min
Sample Name :	Flagelin22			B.conc	:50%H2O/50%MeOH
Cat No. :RP19986					
Lot No. : E1810925H					
Theoretical MW : 2272.50					
Observed MW : 2272.2					