

Sample Name:KRAS G12D Peptide (VVVGADGVGK)
Cat No. : RP30227
Lot No. : E1835326C
Time Processed :19:57:19
Month-Day-Year Processed :03/26/2026


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ25010348

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.233	6055	1006	0.272
2	7.692	6364	1096	0.286
3	8.950	2150014	327964	96.624
4	9.283	45669	6855	2.052
5	16.492	17032	2399	0.765
Total		2225134	339320	100.000