

Sample Name :KRAS G12D Peptide (VVVGADGVGK)

Cat No. : RP30227

Lot No. : E2268724L

Time Processed :12:30:09

Month-Day-Year Processed :01/02/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ19010142

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.825	28402	1469	1.278
2	8.400	16439	1006	0.740
3	9.158	9362	1193	0.421
4	9.433	2145253	394196	96.516
5	9.675	23247	3947	1.046
Total		2222702	401811	100.000