

Sample Name :C-Peptide (3-33), human

Cat No. :RP11116

Lot No. :E2953525J

Time Processed :8:42:28

Month-Day-Year Processed :10/23/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ17010509

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	13.775	6616	1227	0.249
2	14.167	2534127	307244	95.270
3	14.483	16574	3330	0.623
4	14.767	102613	11237	3.858
Total		2659930	323039	100.000