

Sample Name :C-Peptide (3-33), human
 Cat No. :RP11116
 Lot No. :E3744826C
 Time Processed :0:27:10
 Month-Day-Year Processed :03/22/2026
 Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm
 <<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ25010053

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	14.733	12263	1403	0.518
2	15.058	2251152	251112	95.093
3	15.392	13293	3879	0.562
4	15.617	90601	8689	3.827
Total		2367309	265083	100.000