


Mass Spectrum



Sample Information		Interface		Equipment	
Month-Day Processed :	07/07/26	Nebulizing Gas Flow :	1.5L/min	Interface Bias	:ZJ22010150
Time Processed :	9:53:15 AM	CDL Temp :	250	Drying Gas Flow	: +4.5 kV
Injection Volume :	0.4	Block Temp :	200	T.Flow	:5 L/min
Sample Name :	Flagelin22			B.conc	:0.2 ml/min
Cat No. :	RP19986				:50%H2O/50%MeOH
Lot No. :	E4115426G				
Theoretical MW:	2272.50				
Observed MW :	2272.5				