

Sample Name:KRAS G12D Peptide (VVVGADGVGK)
Cat No. : RP30227
Lot No. : E4350625L
Time Processed :2:38:04
Month-Day-Year Processed :12/13/2025


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-3 4.6 x 250 mm
Equipment: SS-CM-0310

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.783	30507	3678	0.660
2	11.200	4396672	730752	95.140
3	11.375	111656	20526	2.416
4	11.583	47059	7572	1.018
5	18.600	35383	5866	0.766
Total		4621278	768394	100.000