

Sample Name: C-Peptide (3-33), human

Cat No. : RP11116

Lot No. : E4633425I

Time Processed :17:16:55

Month-Day-Year Processed :09/25/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ21010375

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	9.250	53919	7557	0.740
2	12.200	25130	3675	0.345
3	12.633	6992022	1118889	95.977
4	13.167	23415	3200	0.321
5	13.550	30430	4274	0.418
6	13.725	160199	21906	2.199
Total		7285116	1159500	100.000