

Sample Name :KRAS G12D Peptide (VVVGADGVGK)

Cat No. : RP30227

Lot No. : E5126226C

Time Processed :19:27:08

Month-Day-Year Processed :04/01/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ20010264

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	5.217	14002	2725	1.031
2	7.933	16271	2723	1.199
3	8.275	5093	1129	0.375
4	8.750	8661	1544	0.638
5	9.075	1294961	218585	95.394
6	9.450	9217	1631	0.679
7	9.667	9276	1522	0.683
Total		1357480	229860	100.000