

Sample Name: C-Peptide (3-33), human
Cat No. :RP11116
Lot No. :E6526625L
Time Processed: 14:23:47
Month-Day-Year Processed :12/07/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm
<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column : Inertsil ODS-SP 4.6 x 250 mm

Equipment:ZJ19010014

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	5.456	1563	176	0.028
2	11.157	33653	2964	0.598
3	11.851	5259	849	0.093
4	17.673	5469871	679312	97.176
5	17.896	49442	8315	0.878
6	18.268	7053	910	0.125
7	18.659	43634	5321	0.775
8	19.052	18331	1811	0.326
Total		5628806	699657	100.000