


Mass Spectrum



Sample Information

Month-Day Processed : 01/30/26
 Time Processed : 9:20:57 AM
 Injection Volume : 0.4
 Sample Name : Humanin, HN
 Cat No. : RP30893
 Lot No. : E7749826A

Theoretical MW: 2687.27

Observed MW : 2686.8

Interface

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment

Equipment : ZI21010035
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50%H2O/50%MeOH