

Sample Name :KRAS G12D Peptide (VVVGADGVGK)
Cat No. : RP30227
Lot No. : E7912925G
Time Processed :2:30:53
Month-Day-Year Processed :07/24/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ17010506

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.083	3553623	382213	96.931
2	6.433	14910	2271	0.407
3	6.658	47599	6364	1.298
4	7.108	28335	3434	0.773
5	11.042	21683	2961	0.591
Total		3666150	397242	100.000