

Sample Name :KRAS G12D Peptide (VVVGADGVGK)

Cat No. : RP30227

Lot No. : E7929525G

Time Processed :19:24:43

Month-Day-Year Processed :07/25/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ25010053

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.017	5776862	526990	97.786
2	6.525	75569	4776	1.279
3	6.933	20452	3303	0.346
4	7.108	34797	5109	0.589
Total		5907681	540178	100.000