

Sample Name : KRAS G12D Peptide (VVVGADGVGK)

Cat No. : RP30227

Lot No. : E8431623E

Time Processed :14:49:45

Month-Day-Year Processed :05/28/2023

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: GR11010440

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	9.902	39339	3655	0.297
2	12.255	3099	316	0.023
3	12.991	45150	5469	0.341
4	13.148	101243	13225	0.765
5	13.393	69780	1390	0.527
6	13.956	12800738	1495161	96.725
7	14.317	128695	6890	0.972
8	14.793	14946	1083	0.113
9	15.965	19057	2149	0.144
10	16.488	8722	572	0.066
11	17.341	3378	94	0.026
Total		13234148	1530006	100.000