

Version 1.0

Update: 07/15/2019

MSH6 Antibody, mAb, Rabbit

Cat. No.	Name	Clone
V08201	MonoRab™ MSH6 Antibody (3B6), mAb, Rabbit	3B6
Specificity	Human MSH6	
Isotype	Rabbit IgG for mAb 3B6	
Production	Cultured <i>in vitro</i> under conditions free from animal-derived components	
Purification	Protein A/G affinity column	
Formulation	50 mM Na-citrate, 150 mM NaCl, pH 7.0, containing 0.03% Proclin 300 *	
Storage	For long term storage, aliquot and store at -20°C or below. Avoid repeated freezing and thawing cycles.	
Application	Platform	Clone
	IHC	3B6
Background	MSH6, also known as MutS protein homolog 6, is a member of the Mutator S (MutS) family of proteins. It functions to repair mutations during DNA replication. MSH6 is a useful in diagnostics of gastrointestinal tract cancer.	
Note	GenScript can customize this product per customer's request including product size, buffer components, etc.	

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.