


# Mass Spectrum

Positive  
Intensity



## Sample Information

User : Gary  
Month-Day Processed : 01/16/14  
Time Processed : 13:18:38  
Injection Volume : 0.2  
Sample Name :  $\beta$ -Amyloid (1-42), human  
Cat No. : RP10017  
Lot No. : P16501407  
Theoretical MW : 4514.10  
Observed MW : 4514.75

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH