


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 01/28/14
Time Processed : 18:37:08
Injection Volume : 0.2
Sample Name : β -Amyloid (1-40)
Cat. No.: RP10004
Lot. No.: P17661405
Theoretical MW : 4329.82
Observed MW : 4330.25

Interface :ESI
Nebulizing Gas Flow :1.5L/min
CDL Temp :250°C
Block Temp :200°C

Interface Bias : +4.5 kV
Drying Gas Flow :5 L/min
T.Flow :0.2 ml/min
B.conc :50%H₂O/50%MeOH