

**Abz-FR-K (Dnp)-P-OH****Cat. No.:** RP20115**Size:** 1 mg**Description:**

Excellent angiotensin I-converting enzyme (ACE) substrate with a Km value of 4.0  $\mu\text{M}$  and a kcat value of 210s<sup>-1</sup>.

**N-Terminal:** Abz**Sequence (one-letter code):**

Abz-FR-K (Dnp)-P-OH

**Sequence (three-letter code):**

Abz-Phe-Arg-Lys(Dnp)-Pro-OH

**Formula:** C<sub>39</sub>H<sub>49</sub>N<sub>11</sub>O<sub>10</sub>**Molecular Weight:** 831.89**Purity:** 95%**Storage:**

Store GenScript Abz-FR-K (Dnp)-P-OH at -20°C.

**Note:**

\*For Non-Clinical Research Use Only \*

