

Rev02  
Update: Nov,23,2023

**DATASHEET**

# Cys-TAT(47-57)

Cat. No.: RP20343

## Overview

<b>Description</b>	A TAT(Transactivated-transcripton) peptide derived from the HIV TAT protein has the ability to transduce peptides or proteins into various cells.
<b>Sequence</b>	{CYS}{TYR}{GLY}{ARG}{LYS}{LYS}{ARG}{ARG}{GLN}{ARG}{ARG}{ARG}
<b>Sequence Shortening</b>	CYGRKKRRQRRR
<b>C Terminal</b>	NH <sub>2</sub>
<b>Molecular Weight</b>	1661.99

## Properties

<b>Purity</b>	> 95%
<b>Form</b>	Lyophilized
<b>Storage</b>	Store at -20°C.

生产商: 南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路28号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China