

Rev01
Update: Sep,23,2022

DATASHEET

Tripeptide Library Matrix (8000 peptide)

Cat. No.: RP30222

Overview

Description

This peptide library is a combination of 8000 (20*20*20) tripeptides. Each tripeptide consists of three amino acid residues out of 20 nature amino acids. This product is delivered in 80 racks of individual 100 peptides in each. GenScript's tripeptide library can be used for different downstream applications as below:

- Ligand-receptor-interactions
- Epitope-mapping
- Vaccine-development
- Drug-screening
- Conformational-epitope analysis
- Spacer-amino-acid positions
- Food-production

Properties

Length	3 aa
Purity	>90%
Form	Lyophilized
Storage	Store at -20°C.
Note	The peptides of this product are supplied as trifluoroacetate salts.