

Rev02  
Update: Aug,24,2023

**DATASHEET**

# SARS-CoV-2 Spike-derived Peptide (YLQPRTFLL)

Cat. No.: RP30230

## Overview

<b>Description</b>	Antibody screening, T-cell assays, Immune monitoring, Antigen specific T-cell stimulation, Cellular immune response
<b>Sequence</b>	{Tyr}{Leu}{Gln}{Pro}{Arg}{Thr}{Phe}{Leu}{Leu}
<b>Sequence Shortening</b>	YLQPRTFLL
<b>Molecular Weight</b>	1150.38

## Properties

<b>Purity</b>	>95%
<b>Solubility</b>	Can be dissolved in DMSO
<b>Form</b>	Lyophilized
<b>Storage</b>	Store at -20°C. Keep tightly closed. Store in a cool dry place.

**For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.**

生产商：南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路28号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China