

## CERTIFICATE OF ANALYSIS

### Product Information

---

Product Name:	Angiostatin K1-3, Human
Cat. No.:	Z02730
Lot No.:	R50451204
Size:	50µg/vial
Storage Condition:	Upon reconstitution, store at 2-8°C for up to one week, and store at -20°C to -80°C for long term storage. Avoid repeated freeze/thaw cycles.

---

### Quality Control

---

Test Items	Specifications	Results
Physical Appearance	White Freeze-dried Powder	Conformed
Purity (SDS-PAGE)	> 95%	Conformed
Purity (HPLC)	> 95%	95.81%
Specific Activity	> $5.5 \times 10^5$ IU/mg	$5.77 \times 10^5$ IU/mg
Molecular Weight	30.0±3.0KDa	Conformed
Bacterial Endotoxin	<0.2EU/µg	<0.2EU/µg

---

### Caution:

For laboratory use only. Not for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:  Bert Date: 04/01/2012  
New Technology Department