

MonoRab™ SARS-CoV-2 Nucleocapsid Antibody (N34), mAb, Rabbit

Cat. No.: A02136

Overview

Specificity	GenScript MonoRab™ SARS-CoV-2 Nucleocapsid Antibody (N34), mAb, Rabbit is specific for SARS-CoV-2 Nucleocapsid protein of wild type, delta variant, and omicron variant.
Host Species	Rabbit
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid protein
Conjugate	Unconjugated

Applications

Working concentrations for specific applications should be determined by the investigators. The appropriate concentrations may be affected by secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.


Application	Recommended Usage
ELISA/CLIA	1:1,000 - 1:5,000
Western Blot	1:100 - 1:1,000
LFA Assay	A02135(C) - A02136 (D)
Other applications	User-optimized

Properties

Form	Liquid
Storage Buffer	PBS, pH 7.2, 0.02% sodium azide
Concentration	1 mg/ml

Storage Instructions	Store at -20°C. This product is stable for 1 year upon receipt, when handled and stored as instructed. Avoid repeated freezing and thawing cycles.
Purification	Affinity chromatography
Isotype	Rabbit IgG
Clonality	Monoclonal
Clone ID	N34
Note	GenScript can customize this product per customer's request including product size, buffer components, etc.


Examples



BLI (Biolayer interferometry) binding affinity measurements of MonoRab™ SARS-CoV-2 Nucleocapsid Antibody (N34), mAb, Rabbit (GenScript, A02136) to nucleocapsid protein of SARS-CoV-2 omicron variant.


MonoRab™ SARS-CoV-2 Nucleocapsid Antibody (N34) captured on Protein A sensors surface can bind omicron variant with dissociation constant (KD) of 6.73E-10M.

SARS-CoV-2 nucleocapsid protein (wild-type)




BLI (Biolayer interferometry) binding affinity measurements of MonoRab™ SARS-CoV-2 Nucleocapsid Antibody (N34), mAb, Rabbit (GenScript, A02136) to SARS-CoV-2 Nucleocapsid protein of wild type and delta variant.

SARS-CoV-2 nucleocapsid protein (delta variant)



MonoRab™ SARS-CoV-2 Nucleocapsid Antibody (N34) captured on Protein A sensors surface can bind wild-type or delta variant with dissociation constant (KD) of 2.25E-09M and 6.83E-10M.

LFA Analysis with A02135 (N338) as capture antibody and A02136 (N34) as detection antibody




Standard curve of SARS-CoV-2 Nucleocapsid Protein LFA Assay.

The LFA assay is developed by using SARS-CoV-2 Nucleocapsid Antibody (N338), mAb, mouse(GenScript, A02135) and MonoRab™ SARS-CoV-2 Nucleocapsid Antibody (N34), mAb, Rabbit(GenScript, A02136) as the capture and detection antibodies, respectively.

The combination shows good dose response of fluorescence intensity and concentrations of SARS-CoV-2 Nucleocapsid protein in LFA assay.

Nucleocapsid protein concentration was nicely correlated with

fluorescence intensity.



Western Blot of wild-type and delta variant Nucleocapsid protein cell lysates with MonoRab™ SARS-CoV-2 Nucleocapsid Antibody (N34), mAb, Rabbit (GenScript, A02136).

M: Protein Marker (GenScript, M00521)

Lane 1: 3 µg Wild-type N protein cell lysate

Lane 2: 3 µg Delta N protein cell lysate

Primary Antibody:

MonoRab™ SARS-CoV-2 Nucleocapsid Antibody (N34), mAb, Rabbit (GenScript, A02136) 1 µg/ml

Secondary Antibody:

Goat Anti-Rabbit IgG-HRP (GenScript, Cat No. A00131)

Background

Target Background : SARS-CoV-2 (Severe acute respiratory syndrome coronavirus 2), also known as 2019-nCoV, is a positive-sense single-stranded RNA virus. It caused coronavirus disease 2019 (COVID-19). Nucleocapsid Protein is a most abundant structure protein of the coronavirus which is associated with nucleic acid.

Synonyms : Rabbit Anti-COVID19 N Protein Antibody; Rabbit Anti-SARS-CoV-2 NP Antibody; Sars/sars-CoV-2 Coronavirus Nucleocapsid Monoclonal Antibody

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.