

Rev03  
 Update: Aug,08,2025
**DATASHEET**


# CDH17/Cadherin 17 Domain 3&4 mFc Chimera, Human

Cat. No.: Z05149


## Product Introduction

<b>Species</b>	Human
<b>Protein Construction</b>	<div style="display: flex; justify-content: space-between; align-items: center;"> <div style="background-color: #0056b3; color: white; padding: 5px; text-align: center;"> <b>CDH17/Cadherin 17 Domain 3&amp;4 (Val245-Phe449)</b>            Accession # Q12864         </div> <div style="background-color: #76b82a; color: white; padding: 5px; text-align: center;"> <b>mFc (IgG1)</b> </div> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> <span>N-term</span> <span>C-term</span> </div>
<b>Purity</b>	> 95% as determined by BisTris PAGE > 95% as determined by HPLC
<b>Endotoxin Level</b>	Less than 1EU per µg by the LAL method.
<b>Expression System</b>	HEK293
<b>Theoretical Molecular Weight</b>	49.16 kDa
<b>Apparent Molecular Weight</b>	Due to glycosylation, the protein migrates to 55-70 kDa based on Bis-Tris PAGE result.
<b>Formulation</b>	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4).
<b>Reconstitution</b>	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Storage &amp; Stability</b>	Upon receiving, the product remains stable up to 6 months at -20 °C or below. Upon reconstitution, the product should be stable for 3 months at -80 °C. Avoid repeated freeze-thaw cycles.

## Examples



CDH17/Cadherin 17 Domain 3&4 mFc Chimera, Human on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.



The purity of CDH17/Cadherin 17 Domain 3&4 mFc Chimera, Human is greater than 95% as determined by SEC-HPLC.

## Background

**Target Background :** Liver-intestine cadherin (CDH17) has been known to function as a tumor stimulator and diagnostic marker for almost two decades. In vivo studies showed CDH17 knockout resulted in apoptotic PC tumor death through activating caspase-3 activity. Taken together, CDH17 functions as an oncogenic molecule critical to PC growth by regulating tumor apoptosis signaling pathways and CDH17 could be targeted to develop an anti-PC therapeutic approach.

**Synonyms :** Cadherin-17; BILL-cadherin; LI-cadherin; P130; Cdh17; CDH16; HPT-1; cadherin-16;

**For research use only. Not intended for human or animal clinical trials, therapeutic or diagnostic use.**

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28 Yongxi Road, Jiangning District, Nanjing, Jiangsu, China