

Rev03
 Update: Aug,08,2025

DATASHEET


LAMP5, His, Human

Cat. No.: Z05560


Product Introduction

Species	Human
Protein Construction	<div style="display: flex; align-items: center; justify-content: center;"> <div style="background-color: #0056b3; color: white; padding: 5px; text-align: center;"> LAMP5 (Glu30-Glu235)_x000D_ Accession # Q9UJQ1-1 </div> <div style="background-color: #76b82a; color: white; padding: 5px; text-align: center; margin-left: 10px;"> His </div> </div> <div style="display: flex; justify-content: space-around; margin-top: 5px; font-size: small;"> N-term C-term </div>
Purity	> 95% as determined by BisTris PAGE > 95% as determined by HPLC
Endotoxin Level	Less than 1EU per µg by the LAL method.
Expression System	HEK293
Theoretical Molecular Weight	24.17 kDa
Apparent Molecular Weight	Due to glycosylation, the protein migrates to 30-40 kDa based on Bis-Tris PAGE result.
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4).
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage & Stability	Upon receiving, the product remains stable up to 6 months at -20 °C or below. Upon reconstitution, the product should be stable for 3 months at -80 °C. Avoid repeated freeze-thaw cycles.

Examples



LAMP5, His, Human on Bis-Tris PAGE under reduced condition.
 The purity is greater than 95%.



The purity of LAMP5, His, Human is greater than 95% as determined by SEC-HPLC.

Background

Target Background : Lysosome-associated membrane protein 5 (LAMP5) is a mammalian ortholog of the *Caenorhabditis elegans* protein, UNC-46, which functions as a sorting factor to localize the vesicular GABA transporter UNC-47 to synaptic vesicles. LAMP5 deficiency led to a larger intensity-dependent increase of wave I, II and V peak amplitude of auditory brainstem response. LAMP5 plays a pivotal role in sensorimotor processing in the brainstem and spinal cord.

Synonyms : BAD-LAMP; LAMP-5; BADLAMP; C20orf103; dJ1119D9.3; UNC-43

For research use only. Not intended for human or animal clinical trials, therapeutic or diagnostic use.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China