

Rev03  
 Update: Aug,08,2025

**DATASHEET**

# CD3G/CD3 gamma, His, Human

Cat. No.: Z04008


## Product Introduction

<b>Species</b>	Human
<b>Protein Construction</b>	<div style="display: flex; align-items: center; justify-content: center;"> <div style="background-color: #0056b3; color: white; padding: 5px; text-align: center;"> <b>CD3G/CD3 gamma (Gln23-Ser116)</b>            Accession # P09693         </div> <div style="background-color: #76923c; color: white; padding: 5px; text-align: center; margin-left: 10px;"> <b>His</b> </div> </div> <div style="display: flex; justify-content: space-around; margin-top: 5px;"> <span>N-term</span> <span>C-term</span> </div>
<b>Purity</b>	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC
<b>Endotoxin Level</b>	Less than 1EU per µg by the LAL method.
<b>Expression System</b>	HEK293
<b>Theoretical Molecular Weight</b>	12.31 kDa
<b>Apparent Molecular Weight</b>	Due to glycosylation, the protein migrates to 25-35 kDa based on Bis-Tris PAGE result.
<b>Formulation</b>	Lyophilized from a 0.22 µm filtered solution in PBS, (pH 7.4).
<b>Reconstitution</b>	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Storage &amp; Stability</b>	Upon receiving, the product remains stable up to 6 months at -20 °C or below. Upon reconstitution, the product should be stable for 3 months at -80 °C. Avoid repeated freeze-thaw cycles.

## Examples



CD3G/CD3 gamma, His, Human on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.



The purity of CD3G/CD3 gamma, His, Human is greater than 95% as determined by SEC-HPLC

## Background

**Target Background :** CD3 gamma, a subunit of the T cell receptor-CD3 (TCR/CD3) complex, helps to support surface TCR/CD3 expression and participates in signal transduction for gene induction after antigen recognition by T lymphocytes, and in TCR/CD3 down-modulation.

**Synonyms :** CD3g; T3G; CD3 gamma; IMD17

**For research use only. Not intended for human or animal clinical trials, therapeutic or diagnostic use.**

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China