

Rev03
Update: Aug,08,2025**DATASHEET**


Siglec-2/CD22, His, Human

Cat. No.: Z05821


Product Introduction

Species	Human
Protein Construction	Siglec-2/CD22 (Asp20-Arg687) Accession # P20273-1 His N-termC-term
Purity	> 95% as determined by Bis-Tris PAGE
Endotoxin Level	Less than 1EU per µg by the LAL method.
Biological Activity	Measured by its binding ability in a functional ELISA. Immobilized Siglec-2/CD22, His, Human at 1µg/ml (100µl/well) on the plate can bind Anti-Siglec-2 Antibody, hFc Tag. Test result was comparable to standard batch.
Expression System	HEK293
Theoretical Molecular Weight	76.2 kDa
Apparent Molecular Weight	Due to glycosylation, the protein migrates to 100-120 kDa based on Bis-Tris PAGE result.
Formulation	Lyophilized from a 0.22 µm filtered solution in PBS, (pH 7.4).
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage & Stability	Upon receiving, the product remains stable up to 6 months at -20 °C or below. Upon reconstitution, the product should be stable for 3 months at -80 °C. Avoid repeated freeze-thaw cycles.

Examples



Immobilized Siglec-2/CD22, His, Human, His Tag at 1 µg/ml (100 µL/well) on the plate. Dose response curve for Anti-Siglec-2 Antibody, hFc Tag with the EC₅₀ of 6.6 ng/ml determined by ELISA.



Siglec-2/CD22, His, Human on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

Background

Target Background : CD22, or cluster of differentiation-22, is a molecule belonging to the SIGLEC family of lectins. It is found on the surface of mature B cells and to a lesser extent on some immature B cells. CD22 a member of the immunoglobulin superfamily. CD22 functions as an inhibitory receptor for B cell receptor (BCR) signaling. It is also involved in the B cell trafficking to Peyer's patches in mice.

Synonyms : CD22 molecule; CD22; Siglec-2; B-cell receptor CD22; BL-CAM; SIGLEC2; SIGLEC2FLJ22814

For research use only. Not intended for human or animal clinical trials, therapeutic or diagnostic use.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China