

Rev03
 Update: Aug,08,2025

DATASHEET


HLA-A*03:01&B2M&KRAS WT (VVVGAGGVGK) Monomer, His & Avi, Human

Cat. No.: Z06607


Product Introduction

Species	Human
Protein Construction	<div style="display: flex; align-items: center;"> <div style="background-color: #0056b3; color: white; padding: 5px; margin-right: 10px;"> HLA-A*03:01&B2M&KRAS WT (VVVGAGGVGK) Monomer[Gly25-Thr305(HLA-A*03:01), Ile21-Met119(B2M) and VVVGAGGVGK peptide] Accession # NP_002107.3(HLA-A*03:01)&P61769(B2M)&VVVGAGGVGK </div> <div style="display: flex; gap: 5px;"> <div style="background-color: #90c090; padding: 5px;">His</div> <div style="background-color: #4f81bd; padding: 5px;">Avi</div> </div> </div> <p style="margin-top: 5px; font-size: small;">N-term C-term</p>
Purity	> 95% as determined by BisTris PAGE > 95% as determined by HPLC
Endotoxin Level	Less than 1EU per µg by the LAL method.
Biological Activity	Measured by its binding ability in a functional ELISA. Immobilized HLA-A*03:01&B2M&KRAS WT (VVVGAGGVGK) Monomer, His & Avi, Human at 0.5µg/ml (100µl/well) on the plate can bind AntiB2M Antibody, mFc Tag. Test result was comparable to standard batch.
Expression System	HEK293
Theoretical Molecular Weight	50.2 kDa
Apparent Molecular Weight	Due to glycosylation, the protein migrates to 52-62 kDa based on Bis-Tris PAGE result.
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4).
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage & Stability	Upon receiving, the product remains stable up to 6 months at -20 °C or below. Upon reconstitution, the product should be stable for 3 months at -80 °C. Avoid repeated freeze-thaw cycles.


Examples



HLA-A*03:01&B2M&KRAS WT (VVGAGGVGK) Monomer, His & Avi, Human on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.



Immobilized HLA-A*03:01&B2M&KRAS WT (VVGAGGVGK) Monomer, His & Avi, Human, His Tag at 0.5 µg/ml (100 µl/well) on the plate. Dose response curve for Anti-B2M Antibody, mFc Tag with the EC50 of 11.2 ng/ml determined by ELISA.



The purity of HLA-A*03:01&B2M&KRAS WT (VVGAGGVGK) Monomer, His & Avi, Human is greater than 95% as determined by SEC-HPLC.

Background

Target Background : Kirsten rat sarcoma 2 viral oncogene homolog (KRAS) is the most commonly mutated oncogene in human cancer. The developments of many cancers depend on sustained expression and signaling of KRAS, which makes KRAS a high-priority therapeutic target.

Synonyms : MHC; KRAS; K-Ras 2; KRAS2; C-K-RAS; CFC2; K-RAS2A; K-RAS2B; K-RAS4A; K-RAS4B; KRAS1; KRAS2; NS; NS3; RASK2; GTPase Kras; KI-RAS; RALD

For research use only. Not intended for human or animal clinical trials, therapeutic or diagnostic use.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China