

Rev03
Update: Aug,08,2025**DATASHEET**

MonoRab™ Anti-Lecanemab Antibody (44G1), mAb, Rabbit

Cat. No.: A02270

Overview

Specificity	This product is specific for Lecanemab.
Host Species	Rabbit
Immunogen	Lecanemab
Conjugate	Unconjugated

Applications

Working concentrations for specific applications should be determined by the investigators. The appropriate concentrations may be affected by secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.

Application	Recommended Usage
Sandwich ELISA	0.5-2 µg/ml
ELISA	0.01-1 µg/ml
Competitive ELISA	1-10 µg/ml


Properties

Form	Lyophilized
Storage Buffer	Lyophilized with PBS, pH 7.2, containing 0.02% sodium azide.
Reconstitution	Reconstitute the lyophilized antibody with deionized water (or equivalent) to a final concentration of 0.5 mg/mL.
Storage Instructions	The lyophilized product remains stable for up to 1 year at -20 °C from the date of receipt. Upon reconstitution, it can be stored for 2-3 weeks at 2-8 °C or for up to 12 months at -20 °C or below. Avoid repeated freezing and thawing cycles.
Purification	Protein A affinity column

Isotype	Rabbit IgG, κ
Clonality	Monoclonal
Clone ID	44G1
Note	GenScript can customize this product per customer's request including product size, buffer components, etc.

Examples


Anti-Lecanemab Antibody (44G1) or (4F2) inhibits the binding of Lecanemab to β -Amyloid (1-42)



Anti-Lecanemab Antibody (44G1) or (4F2) inhibits the binding of Lecanemab to β -Amyloid (1-42). Coating antigen: β -Amyloid (1-42), human (GenScript, RP10017) 1 μ g/ml.

Lecanemab final concentration: 10 μ g/ml.


Anti-Lecanemab antibody dilutions start from 10 μ g/ml. MonoRab™ Anti-Lecanemab Antibody (44G1), mAb, Rabbit (GenScript, A02270) and MonoRab™ Anti-Lecanemab Antibody (4F2), mAb, Rabbit (GenScript, A02271) demonstrated inhibitory properties.



BLI (Biolayer interferometry) binding affinity measurements of MonoRab™ Anti-Lecanemab Antibody (44G1), mAb, Rabbit (GenScript, A02270) to Lecanemab.

Lecanemab captured on HFC (Anti-Human IgG Fc) Probes can bind MonoRab™ Anti-Lecanemab Antibody (44G1), mAb, Rabbit with a dissociation constant (KD) of less than 1 pM.

MRD Analysis of the Detection of Lecanemab in Biological Matrix




MRD Analysis of the Detection of Lecanemab in Biological Matrix.

The MRD is the minimum dilution necessary for the detection of Lecanemab in biological matrix with least interference. Serum samples from cynomolgus monkey were serially diluted to determine the MRD of this assay, and the test result suggested that MRD was as 1:10.

In this ELISA assay, MonoRab™ Anti-Lecanemab Antibody (44G1), mAb, Rabbit (GenScript, A02270) was coated at a concentration of 1 μ g/ml, and MonoRab™ Anti-Lecanemab Antibody (4F2), mAb, Rabbit (GenScript, A02271) conjugated with Biotin was used as a detection antibody at a concentration of 0.5 μ g/mL.

Specificity analysis of the detection of Lecanemab




Specificity analysis of the detection of Lecanemab.

Lecanemab samples at 8 concentrations were spiked with 200 ng/mL of human IgG1. The test result demonstrated that the high concentration of human IgG1 did not interfere with the detection of Lecanemab.

In this ELISA assay, MonoRab™ Anti-Lecanemab Antibody (44G1), mAb, Rabbit (GenScript, A02270) was coated at a concentration of 1 µg/ml, and MonoRab™ Anti-Lecanemab Antibody (4F2), mAb, Rabbit (GenScript, A02271) conjugated with Biotin was used as a detection antibody at a concentration of 0.5 µg/mL.

Standard Curve of Lecanemab Sandwich ELISA by A02270 (44G1) and A02271 (4F2)



Standard Curve of Lecanemab Sandwich ELISA by A02270 (44G1) and A02271 (4F2)

In this Sandwich ELISA assay, MonoRab™ Anti-Lecanemab Antibody (44G1), mAb, Rabbit (GenScript, A02270) was coated at a concentration of 1 µg/ml, and MonoRab™ Anti-Lecanemab Antibody (4F2), mAb, Rabbit (GenScript, A02271) conjugated with Biotin was used as a detection antibody at a concentration of 0.5 µg/mL. GenScript can provide customized conjugation service for this product per customer's request. In this assay, a four-parameter logistic curve fitting program was used to create a standard curve with the R-Square equal to 0.99971. The typical dynamic range of the assay is 0.31 - 20 ng/mL and its detection limit is 0.31 ng/mL.

Background

Target Background : Lecanemab, marketed as Leqembi, is a medication that is prescribed to treat Alzheimer's disease. Lecanemab is a monoclonal antibody derived from a mouse antibody called mAb158, which has been humanized to improve its compatibility with humans. It specifically targets amyloid beta, a protein associated with the disease.

Synonyms : Lecanemab; Leqembi; Lecanemab-irmb; BAN2401; mAb158; Anti-amyloid beta antibody

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China