


Rev03
 Update: Aug,08,2025

DATASHEET


NKG2C/CD159c[Biotin], His & Avi, Human

Cat. No.: Z04569

Product Introduction

Species	Human
Protein Construction	
Conjugate	Biotin
Purity	> 95% as determined by BisTris PAGE > 95% as determined by HPLC
Endotoxin Level	Less than 1EU per µg by the LAL method.
Expression System	HEK293
Theoretical Molecular Weight	15.3 kDa
Apparent Molecular Weight	Due to glycosylation, the protein migrates to 35-55 kDa based on Bis-Tris PAGE result.
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4).
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage & Stability	Upon receiving, the product remains stable up to 6 months at -20 °C or below. Upon reconstitution, the product should be stable for 3 months at -80 °C. Avoid repeated freeze-thaw cycles.

Examples



NKG2C/CD159c[Biotin], His & Avi, Human on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

Background

Target Background : As a first line of defense, natural killer (NK) cells play a crucial role in the fight against infections. The presented study is the first of its kind that ascribes CD56dimCD16 NKG2C-expressing NK cells a crucial role in biasing adaptive immune responses upon influenza vaccination and suggests NKG2C as a potential biomarker in predicting pandemic influenza vaccine responsiveness.

Synonyms : CD159c; KLRC2; NKG2C; NK cell receptor C

For research use only. Not intended for human or animal clinical trials, therapeutic or diagnostic use.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28 Yongxi Road, Jiangning District, Nanjing, Jiangsu, China