


Rev07  
Update: Aug,11,2025**DATASHEET**

# CD24 Fc Chimera, Human

Cat. No.: Z03568

## Product Introduction

<b>Species</b>	Human
<b>Protein Construction</b>	 CD24 (Ser27-Gly59) Accession # P25063      hFc N-term                                  C-term
<b>Purity</b>	> 90% as analyzed by SDS-PAGE
<b>Endotoxin Level</b>	< 1 EU/μg of protein by gel clotting method
<b>Expression System</b>	HEK 293
<b>Apparent Molecular Weight</b>	~42 kDa, on SDS-PAGE under reducing conditions.
<b>Formulation</b>	Lyophilized from a 0.2 μm filtered solution in PBS.
<b>Reconstitution</b>	It is recommended that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute the lyophilized powder in ddH <sub>2</sub> O or PBS up to 100 μg/ml.
<b>Storage &amp; Stability</b>	Upon receiving, this product remains stable for up to 6 months at lower than -70°C. Upon reconstitution, the product should be stable for up to 1 week at 4°C or up to 3 months at -20°C. For long term storage it is recommended that a carrier protein (example 0.1% BSA) be added. Avoid repeated freeze-thaw cycles.

## Background

**Target Background :** Signal transducer CD24 also known as cluster of differentiation 24 or heat stable antigen CD24 (HSA) is a protein that in humans is encoded by the CD24 gene. CD24 is a cell adhesion molecule. CD24 is a sialoglycoprotein expressed at the surface of most B lymphocytes and differentiating neuroblasts. It is also expressed on neutrophils and neutrophil precursors from the myelocyte stage onwards. The encoded protein is anchored via a glycosyl phosphatidylinositol (GPI) link to the cell surface. The protein also contributes to a wide range of downstream signaling networks and is crucial for neural development. Cross-linking of CD24 on the surface of neutrophils induces apoptosis, and this appears to be defective in sepsis.

**Synonyms :** CD24A; FLJ22950; FLJ43543; MGC75043

**For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.**

Confidential and Privileged



Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China

---

**GenScript USA, Inc.**

860 Centennial Ave. Piscataway, NJ 08854

Tel: 1-732-885-9188