

Rev07
Update: Aug,11,2025**DATASHEET**

MonoRab™ Anti-Bevacizumab Antibody (30E1), mAb, Rabbit

Cat. No.: A01895

Overview

Specificity	The product is specific for Bevacizumab but can't block vascular endothelial growth factor 165 (VEGF 165). The antibody is recommended as a capture antibody in a pharmacokinetic (PK) bridging assay with detection antibody GenScript, A01896-40, MonoRab™ Anti- Bevacizumab Antibody (46E3)[Biotin], mAb, Rabbit.
Host Species	Rabbit
Immunogen	Bevacizumab
Conjugate	Unconjugated

Applications

Working concentrations for specific applications should be determined by the investigators. The appropriate concentration may be affected by secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.


Application	Recommended Usage
ELISA	0.01-0.1 µg/ml

Properties

Form	Lyophilized
Storage Buffer	lyophilized with PBS, pH 7.4, containing 0.02% sodium azide.
Reconstitution	Reconstitute the lyophilized powder with deionized water (or equivalent) to an final concentration of 0.5 mg/mL.
Storage Instructions	The lyophilized product remains stable up to 1 year at -20 °C from date of receipt. Upon reconstitution, it can be stored for 2-3 weeks at 2-8 °C or for up to 12 months at -20 °C or below. Avoid repeated freezing and thawing cycles.
Purification	Protein A affinity column

Isotype	Rabbit IgG
Clonality	Monoclonal
Clone ID	30E1
Note	GenScript can customize this product per customer's request including product size, buffer components, etc.

Examples




MonoRab™ Anti-Bevacizumab Antibody (30E1), mAb, Rabbit (GenScript, A01895-40) binds with Bevacizumab. While the antibody does not recognize the human IgG Fc fragment (data not shown).

Coating antigen: Bevacizumab, 1 µg/ml.

Anti-Bevacizumab antibody (GenScript, A01895-40) dilution start from 2,000 ng/ml.

EC₅₀= 7.280 ng/ml.



Standard curve of Bevacizumab Sandwich ELISA. The Bevacizumab Sandwich ELISA assay is developed by using MonoRab™ Anti-Bevacizumab Antibody (30E1), mAb, Rabbit (GenScript, A01895-40) and MonoRab™ Anti-Bevacizumab Antibody (46E3) [Biotin], mAb, Rabbit (GenScript, A01896-40) as the capture and detection antibodies, respectively.

The sensitivity of detecting Bevacizumab is up to 20 pg/ml.

Background

Target Background : Bevacizumab (Avastin) is a humanized monoclonal antibody that is approved by the U.S. Food and Drug Administration for the treatment of patients with colon cancer, lung cancer, glioblastoma, renal-cell carcinoma, and age-related macular degeneration. As an angiogenesis inhibitor, Bevacizumab blocks angiogenesis by inhibiting vascular endothelial growth factor A (VEGF-A).

Synonyms : Rabbit Monoclonal to Avastin

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

Confidential and Privileged



Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China

GenScript USA, Inc.

860 Centennial Ave. Piscataway, NJ 08854

Tel: 1-732-885-9188