
Human Recombinant G-Protein Coupled Receptor 103 Stable Cell Line**Cat. No. M00321****Version 07292020**

I. Introduction	1
II. Background	1
III. Representative Data	2
IV. Thawing and Subculturing	2
V. References	3
VI. Limited Use License Agreement.....	4

I. Introduction

Catalog Number: M00321

Cell Line Name: CHO-K1/GPR103/G α 15

Gene Synonyms: AQ27; GPR103; MGC149217; SP9155

Expressed Gene: Genbank Accession Number NM_198179; no expressed tags

Host Cell: CHO-K1/G α 15

Culture Properties: Adherent

Quantity: Two vials of frozen cells (>1 \times 10⁶ per vial)

Stability: More than 16 passages

Application: Functional assay for GPR103 receptor (Calcium flux assay)

Freeze Medium: 45% culture medium, 45% FBS (Cat. #10099-141, Gibco), 10% DMSO (Cat. #D2650, Sigma)

Complete Growth Medium: Ham's F-12K (Kaighn's) (Cat. #21127, Life Technologies), 10% FBS

Culture Medium: Ham's F-12K (Kaighn's), 10% FBS, 100 μ g/ml Hygromycin B (Cat. #10687010, Invitrogen), 400 μ g/ml G418 (Cat. #10131-035, Gibco)

Mycoplasma Status: Negative*

Storage: Liquid nitrogen immediately upon receipt.

II. Background

GPR103 is known as an orphan G protein-coupled receptor with reported expression in brain, heart, kidney, adrenal gland, retina, and testis. GPR103 mRNA has been reported to be highly expressed in the superficial layers of the entire spinal cord and a high density of 26RFa binding sites was observed in the superficial layers of the dorsal horn. QRFP binds and activates the human GPR103, as well as mouse GPR103A and GPR103B, with nanomolar affinities in transfected cells. Therefore, the current experiments investigated the effects of QRFP administration in rats and the effects of a high fat diet on prepro-QRFP mRNA and GPR103 receptor mRNA levels. 26RFa and QRFP are endogenous ligands of GPR103.

* The mycoplasma test was performed with MycoAlert™ PLUS Mycoplasma Detection Kit of Lonza.

III. Representative Data

Concentration-dependent stimulation of intracellular calcium mobilization by Human 26RFa in CHO-K1/GPR103/G α 15 cells




Figure 1. Human 26RFa-induced concentration-dependent stimulation of intracellular calcium mobilization in CHO-K1/GPR103/G α 15 cells. The cells were loaded with Calcium-4 prior to being stimulated with agonist human 26RFa. The intracellular calcium change was measured by FLIPR^{TETRA}. The relative fluorescent units (RFU) were normalized and plotted against the log of the cumulative doses of human 26RFa (Mean \pm SD, n = 2). The EC₅₀ of human 26RFa on this cell was 2.0 nM.

Notes:

1. EC₅₀ value is calculated with four parameter logistic equation:

$$Y = \text{Bottom} + \frac{(\text{Top} - \text{Bottom})}{1 + 10^{((\text{LogEC}_{50} - X) * \text{HillSlope})}}$$

X is the logarithm of concentration. Y is the response
 Y is RFU and starts at Bottom and goes to Top with a sigmoid shape.
2. Signal to background Ratio (S/B) = Top/Bottom

IV. Thawing and Subculturing

Thawing: Protocol

1. Remove the vial from liquid nitrogen tank and thaw the cells quickly in a 37°C water-bath.
2. Just before the cells are completely thawed, decontaminate the outside of the vial with 70% ethanol and transfer the cells to a 15 ml centrifuge tube containing 9 ml of complete growth medium.
3. Pellet cells by centrifugation at 200 x g for 5 min, and discard the medium.
4. Resuspend the cells in complete growth medium.
5. Add 10 ml of the cell suspension in a 10 cm dish.
6. Add Hygromycin B and G418 to concentrations of 100 μ g/ml and 400 μ g/ml respectively the

following day.

Subculturing: Protocol

1. Remove and discard culture medium.
2. Wash cells with PBS (pH=7.4) to remove all traces of serum that contains trypsin inhibitor.
3. Add 2.0 ml of 0.05% (w/v) Trypsin- EDTA (GIBCO, Cat No. 25300) solution to a 10 cm dish and observe the cells under an inverted microscope until cell layer is dispersed (usually within 3 to 5 minutes).

Note: To avoid clumping, do not agitate the cells by hitting or shaking the dish while waiting for the cells to detach. Cells that are difficult to detach may be placed at 37°C to facilitate dispersal.

4. Add 6.0 to 8.0 ml of complete growth medium and aspirate cells by gently pipetting, centrifuge the cells 200 x g for 5min, and discard the medium.
5. Resuspend the cells in culture medium and add appropriate aliquots of the cell suspension to new culture vessels.
6. Incubate cultures at 37°C.

Subcultivation Ratio: 1:3 to 1:8

Medium Renewal: Every 2 to 3 days

V. References

1. Baribault H, Danao J, Gupte J, *et al.* (2006) The G-protein-coupled receptor GPR103 regulates bone formation. *Mol. Cell. Biol.* 26(2):709-17.
2. Takayasu S, Sakurai T, Iwasaki S, *et al.* (2006) A neuropeptide ligand of the G protein-coupled receptor GPR103 regulates feeding, behavioral arousal, and blood pressure in mice. *Proc Natl Acad Sci U S A.* 103(19):7438-43
3. Yamamoto T, Wada T, Miyazaki R. (2008) Analgesic effects of intrathecally administered 26RFa, an intrinsic agonist for GPR103, on formalin test and carrageenan test in rats. *Neuroscience.* 157(1):214-22
4. Primeaux SD, Blackmon C, Barnes MJ, *et al.* (2008) Central administration of the RFamide peptides, QRFP-26 and QRFP-43, increases high fat food intake in rats. *Peptides.* 29(11):1994-2000

GenScript USA Inc,
860 Centennial Ave.
Piscataway, NJ 08854
Toll-Free: 1-877-436-7274
Tel: 1-732-885-9188, Fax: 1-732-210-0262
Email: product@genscript.com
Web: <http://www.genscript.com>

For Research Use Only.

Limited Use License Agreement

This is a legal agreement between you (Licensee) and GenScript USA Inc. governing use of GenScript's stable cell line products and protocols provided to licensee. By purchasing and using the stable cell line, the buyer agrees to comply with the following terms and conditions of this label license and recognizes and agrees to such restrictions:

- 1) The products are not transferable and will be used at the site where they were purchased. Transfer to another site owned by buyer will be permitted only upon written request by buyer followed by subsequent written approval by GenScript.
- 2) The purchaser cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party.
- 3) The products sold by GenScript are for laboratory and animal research purposes only. The products are not to be used on humans, for consumption, or for any unlawful uses.

GenScript USA Inc. will not assert against the buyer a claim of infringement of patents owned or controlled by GenScript USA Inc. and claiming this product based upon the manufacture, use or sale of a clinical diagnostic, therapeutic and vaccine, or prophylactic product developed in research by the buyer in which this product or its components has been employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on the use of this product for other purposes, contact Marketing Department, GenScript USA Inc., 860 Centennial Ave. Piscataway, NJ 08854, U.S.A. Phone: 1-732-885-9188. Fax: 1-732-210-0262. Email: marketing@genscript.com.