

Certificate of Analysis

Product ID: OHu22413D Lot. No.: 774762-401/C84565

Product Information:

Gene Symbol: SLC25A4 Accession No.: NM_001151

Cloning Vector: pcDNA3.1⁺/C-(K)DYK Gene Length*: 897 bp

Cloning Strategy: CloneEZ Plasmid size: 6287 bp

Antibiotic resistance: Ampicillin, recommended concentration: 50µg/ml Amount: 10 µg

Note: * The stop codon has not been included in gene length calculation. And the detail sequence showed in attachment 2

Storage: Store at -20°C
The ORF clone is shipped at room temperature, but should be kept at -20°C for long-term storage.

QC Items	Specifications	Results	
Sequencing Alignment	Sequencing results are consistent with the targeted insert sequence.	Pass	Consistent
Vector Sequence	The flanking sequences of the cloning sites are correct.	Pass	Correct Shown in the SQD file
Restriction Digests	The size of inserted fragment is correct and free of unexpected bands suggesting contamination.	Pass	Correct Shown in attachment 1
DNA Quality	Miniprep: 10 µg OD260/280= 1.8~2.0 Free of contamination	Pass	≥ 10 µg OD260/280=1.87 Pure
Quality Grade	Research Grade	Pass	Research Grade
Appearance	Clear and free of foreign particles.	Pass	Clear Free of foreign particles
Additional Test		N/A	

Certified by *Morgan* Date: 12/15/2015

Thank you for your patronage to our Gene services! To maintain this working relationship, we shall be grateful if you can add our webpage URL into your lab website. As a token of appreciation, you will be rewarded by 1,000 EZcoupon™ points. For more information, please contact us by e-mail at web@genscript.com.

For research use only

860 Centennial Ave., Piscataway, NJ 08854, USA

Toll-Free: 1-877-436-7274

Tel: 1-732-885-9188


Fax: 1-732-210-0262

Email: order@genscript.com


Web: www.genscript.com

Plasmid Construct Map

The ORF was cloned in pcDNA3.1⁺/C-(K)DYK by CloneEZ.




Enzyme Digestion



Lane M: KB Ladder
Lane 1: OHu22413D plasmid
Lane 2: OHu22413D plasmid digested by HindIII and SmaI

Digestion Conditions:
~ 300ng plasmid digested
Digestion in water-bath, 37°C for 40 minutes
1% Agarose Gel



For research use only

860 Centennial Ave., Piscataway, NJ 08854, USA

Attachment 2**Original ORF Sequence (Including start and stop codons)**

```
ATGGGTGATC ACGCTTGGAG CTCCTAAAG GACTTCCTGG CCGGGGGCGT CGCCGCTGCC GTCTCCAAGA CCGCGGTCGC CCCCATCGAG AGGGTCAAAC
TGCTGCTGCA GGTCCAGCAT GCCAGCAAAC AGATCAGTGC TGAGAAGCAG TACAAAGGGA TCATTGATTG TGTGGTGAGA ATCCCTAAGG AGCAGGGCTT
CCTCTCCTTC TGGAGGGGTA ACCTGGCCAA CGTGATCCGT TACTTCCCA CCAAGCTCT CAACTTCGCC TTCAAGGACA AGTACAAGCA GCTCTTCTTA
GGGGGTGTGG ATCGGCATAA GCAGTTCTGG CGTACTTTG CTGGTAACCT GCGTCCGGT GGGCCCGCTG GGGCCACCTC CCTTTGCTTT GTCTACCCGC
TGGACTTTC TAGGACCAGG TTGGCTGCTG ATGTGGGCAA GGGCGCCGCC CAGCGTGAGT TCCATGGTCT GGGCGACTGT ATCATCAAGA TCTTCAAGTC
TGATGGCCTG AGGGGGCTCT ACCAGGGTTT CAACGTCTCT GTCCAAGGCA TCATTATCTA TAGAGCTGCC TACTTCGGAG TCTATGATAC TGCCAAGGGG
ATGCTGCCTG ACCCAAGAA CGTGCACATT TTTGTAGCT GGATGATTGC CCAGAGTGTG ACGGCAGTCG CAGGGCTGGT GTCCTACCCC TTTGACTGTG
TTCGTGCTAG AATGATGATG CAGTCCGCC GGAAAGGGC CGATATTATG TACACGGGA CAGTTGACTG CTGGAGGAAG ATTGCAAAG ACGAAGGAGC
CAAGGCCTTC TTCAAAGGTG CCTGGTCAA TGTGCTGAGA GGCATGGCG GTGCTTTGT ATTGGTGTG TATGATGAGA TCAAAAATA TGTCTAA
```

For research use only

860 Centennial Ave., Piscataway, NJ 08854, USA

Toll-Free: 1-877-436-7274

Tel: 1-732-885-9188

Fax: 1-732-210-0262

Email: order@genscript.comWeb: www.genscript.com